


SOLAX MOBILE FOLDABLE BASKET – Model M-FB01-14

Installation Instructions

Package Contents

- A. (1) Foldable Basket
- B. (1) Rear Mounting Bracket Housing
- C. (1) Front Mounting Bracket
- D. (1) Basket Mounting Bracket
- E. (4) Phillips screws




Step 1

Take the basket mounting bracket (D) and place it over the front mounting bracket (C), aligning the four holes. These will get mounted as one unit. Then take the rear mounting bracket housing (B) and attach it to the tiller with the front basket/front mounting bracket as shown.



Step 2

Connect the basket mounting/front mounting bracket with the rear mounting housing with the four Phillips screws (E) provided.



Step 3

Unfold the basket (A) and attached it to the front basket mounting bracket by locking it over the bottom and top mounting handles as shown below. To remove the basket, simply lift the basket up with the handles.



SOLAX
MOBILITY

www.solaxmobility.com

800-983-1306